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GB Directions for use

Photo-activated disinfection of root canals

Composition of **LaserHF ENDO-PDT** solution: Medical grade tolonium chloride powder, water, sodium phosphate buffer addition hydroxy methyl propyl cellulose. 1 ml contains 155 mg tolonium chloride powder.

Indications

This product is for use in disinfection of root canals only.

Contraindications

Known allergic reactions to any of the components of **LaserHF ENDO-PDT** solution.

Side effects

At the concentration provided, **LaserHF ENDO-PDT** solution produces no side effect other than a transient colour change in the dentine on the surface. This is removed by the washing process performed after disinfection of the lesion as described in the application below.

Caution

Contamination with saliva must be avoided to prevent reinfection of the lesion. Avoid contact with eyes but in such cases wash out immediately with water until eye is clear.

Application

Using conventional procedures, isolate tooth using rubber dam and then disinfect the surrounding area and tooth surface with a suitable disinfectant. Access the tooth and remove dentine to expose the pulp chamber. Prepare root canal. Use preferred instrumentation and cleaning agents to remove organic material. Wash out canal thoroughly with sterile solution and ensure no traces of cleansing agents used in preparation remain. Dry canal thoroughly by aspiration followed by sterile paper points. For the application of **LaserHF ENDO-PDT** solution use e.g. Hager & Werken Miraject Endo canula. Insert needle until resistance is felt against canal wall and fill canal and pulp chamber with **LaserHF ENDO-PDT** solution. Agitate the solution in the canal for 60 seconds using an endodontic file one size smaller than that to which the canal has been prepared. For the subsequent laser treatment with PDT/LLP modus following the instructions in the laser operating manual, introduce the fibre into the canal until resistance is encountered. Activate laser according to instructions in Laser Operating Manual for recommended time at recommended power setting while holding fibre in canal. The whole length of the canal is irradiated, to fully activate all **LaserHF ENDO-PDT** solution and to maximise disinfection.

Note: When the distance from the entrance in the pulp chamber to the point of resistance in the canal exceeds 10 mm, treat the canal in sections by withdrawing the fibre in 5 mm steps and then irradiate each section for recommended time. After canal has been completely treated, wash out canal with sterile solution and dry by aspiration followed by sterile paper points. When ready continue with obturation procedure. Fill canal with appropriate sealing agent or dress with medicament according to operator's preference. If a medicament is placed, repeat the **LaserHF ENDO-PDT** procedure at subsequent visit immediately before sealing the canal permanently.

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HAGER & WERKEN GmbH & Co. KG · Ackerstraße 1 · 47269 Duisburg, Germany